Method of Test for DETERMINING ALKALI CONTENT OF FLY ASH BY X-RAY FLUORESCENCE SPECTROMETER DOTD Designation 531

I. Scope

Use this method to determine total sodium oxide (Na₂O) and potassium oxide (K₂O) content of fly ash for calculation of equivalent percentage of sodium oxide by weight.

II. Apparatus

- A. Fly Ash Worksheet (Figure 1)
- B. Analytical Balance 120 g capacity accurate to +/-0.0001 g (Figure 2)
- C. Bead Fusion Instrument (Figure 3)
- D. X-Ray Fluorescence (Figure 4)

III. Reagents

- A. Fusion flux
- B. Release agent

IV. Health Precautions

Proper precautions are to be taken whenever hot materials or equipment must be handled. Use container holder or thermal gloves while handling hot containers. Wear eye protection while stirring and weighing materials. Perform the preparation of sample under vent hood to prevent exposure to fumes.

v. Sample

Obtain sample according to DOTD S102 sampling procedure as stated in Material sampling Manual.

VI. Procedure

- A. Weigh 0.6 g of fly ash and transfer to bead fusion instrument crucible.
- B. Weigh 1.0 g of fusion flux and transfer to bead fusion instrument crucible.
- C. Initiate the fusion procedure in accordance with the instrument's operating instructions.
- D. Upon completion of the fusion cycle, remove the bead from the instrument's mold.
- E. Determine Na₂O and K₂O content utilizing x-ray fluorescence spectrometer in accordance with the instrument manufacturer's operating procedure.

VII. Calculation

- A. Calculate the results as percent by weight of the original fly ash sample.
- B. The equivalent percentage of sodium oxide (Na₂O), calculated as follows:

Equivalent $Na20, \% = Na20, \% + (0.658 \times K20, \%)$

DOTD TR 531 Rev 10/4/2017 Page **2** of **5**

Example:

 $Na_2O = 1.42\%$ $K_2O = 0.17\%$

Equivalent Na₂O % = $1.42 + (0.658 \times 0.17)$

Equivalent Na₂O % = 1.42 + 0.11

Equivalent Na₂O % = 1.53

Equivalent Na₂O % = 1.5

VIII. Report

The data shall be reported to the nearest tenth of a percent (0.1%) on the Fly Ash worksheet (Figure 1).

Example: 1.53 is reported as 1.5

FLY ASH, CLASS F, MATERIAL CODES 297 AND 589

CHEMICAL PROPERTIES

TOTAL OXIDES (ASTM C311)			84.2	Ρ
SILICON DIOXIDE AI2O3	19.9	57.3		
AMMONIUM HYDROXIDE GROUP	7.0	26.9		
SULFUR TRIOXIDE (ASTM C311)		0.3		Ρ
MAGNESIUM OXIDE (ASTM C311)		2.0		
CALCIUM OXIDE (ASTM C311)		2.0		
LOSS ON IGNITION, % (ASTM C311)				
TOTAL ALKALI EQUIVALENT, % (DOTD TR 531)			2.57	
Na2O K2O	1.07 2.28			
MOISTURE CONTENT, % (ASTM C311)				s <u></u>
PASS/FAIL BASED ON AASHTO M295, TABLE 1				
TESTED BY:	I	DATE:	9/28/2017	
CHECKED BY:	I	DATE: _		

085136 F

Figure 1. Fly Ash Worksheet

DOTD TR 531 Rev 10/4/2017 Page **4** of **5**



Figure 2. Analytical Balance



Figure 3. Bead Fusion Instrument

DOTD TR 531 Rev 10/4/2017 Page **5** of **5**



Figure 4. X-Ray Fluorescence